

REICHERT FORESEE PHP™ TECHNICAL SPECIFICATIONS

STIMULUS	Dot deviation signal flashed over macular loci
AREA OF FIELD TESTED	Central visual field - 14°
MODE OF OPERATION	40GB hard drive Removable USB memory stick for results export (formats: HTML, JPEG, BMP) XP embedded software Office network ready
DEVICE DIMENSIONS	Height: 12 in. (30.48 cm) Width: 20 in. (50.8 cm) Depth: 12 in. (30.48 cm)
CHIN REST AND HOME BASE DIMENSIONS	Height: 8in. (20.32 cm) Width: 16in. (40.64 cm) Depth: 16in. (40.64 cm)
OVERALL SYSTEM WEIGHT	40lbs (18.14 kg)
POWER REQUIREMENTS	100-240V, 50/60 Hz, 1A max IEC-320 standard power inlet connector for worldwide use